

For a limited time, 3 for 2 on selected products from the SensiMix range*

Introducing SensiMix™. Simple, Sensitive qPCR Solutions.

Excellent Value

- All-in-one SensiMix kits have been developed for use with Sybr Green I and Probe based chemistry

Excellent Performance

- The SensiMix range has been designed to provide specificity, reproducibility and sensitivity

Excellent Range

- Universal SensiMix range for DNA/cDNA templates
- Premixed SensiMixPlus range for Sybr Green I assays on DNA/cDNA templates
- Universal SensiMix One-Step and Two-Step kits for RNA templates
- SensiMix HRM for high-resolution melt analysis

Please quote LTF/3FOR2 on your order to qualify.

Order now

Please reference Discount Code: **LTF/3FOR2** to validate this offer on orders placed by phone, fax or email.

-lff- Labortechnik GmbH & Co. KG
Hatthauer Str. 18
88142 Wasserburg/B
Deutschland

Tel: +49 (0)8382-9852-0
Fax: +49 (0)8382-9852-32
email: info@labortechnik.com
www.labortechnik.com

Please contact
support@labortechnik.com
for SensiMix and
instrument compatibility.



www.quantace.com

For a limited time, 3 for 2 on selected products from the SensiMix range*

*The offer expires 31st August, 2009. This special offer cannot be used in conjunction with any other offer or discount. Lowest priced item for free. Maximum 5 packs per purchase. To claim your discount, please quote **LTF/3FOR2** on all orders placed by phone, fax or email.

Product Name	Description	Pack Size	Cat. No.	List Price
Universal range of SensiMix reagents for DNA/cDNA templates - suitable for probe or SYBR® Green I assays				
SensiMix dT				
2x SensiMix dT 50mM MgCl ₂ 50x SYBR® Green I	2x formulation, dT based, contains internal reference	250 x 50µl reactions	QT6T3-02	€290.00
		500 x 50µl reactions	QT6T3-05	€550.00
SensiMix dU				
2x SensiMix dU, 50mM MgCl ₂ 50x SYBR® Green I 50x UNG	2x formulation, dU based, contains internal reference	250 x 50µl reactions	QT6U5-02	€320.00
		500 x 50µl reactions	QT6U5-05	€605.00
SensiMix Lite				
5x SensiMix Lite, 50mM MgCl ₂ SensiMix Lite enzyme mix 50x SYBR® Green I	5x formulation, dT based	100 x 20µl reactions	QT405-01	€125.00
		500 x 20µl reactions	QT405-05	€525.00
SensiMix NoRef				
2x SensiMix NoRef 50mM MgCl ₂ 50x SYBR® Green I	2x formulation, dT based, no internal reference	250 x 50µl reactions	QT505-02	€290.00
		500 x 50µl reactions	QT505-05	€550.00
SensiMix LowRef				
2x SensiMix LowRef 50mM MgCl ₂ 50x SYBR® Green I	2x formulation, dT based, low internal reference	250 x 50µl reactions	QT525-02	€290.00
		500 x 50µl reactions	QT525-05	€550.00
All-in-one range of SensiMix reagents for DNA/cDNA templates - exclusively for SYBR® Green I assays				
SensiMixPlus SYBR				
2x SensiMixPlus SYBR 50mM MgCl ₂	2x formulation, dT based, contains internal reference and SYBR® Green I	250 x 50µl reactions	QT605-02	€290.00
		500 x 50µl reactions	QT605-05	€550.00
SensiMixPlus SYBR & Fluorescein				
2x SensiMixPlus SYBR & Fluorescein 50mM MgCl ₂	2x formulation, dT based, contains fluorescein and SYBR® Green I	250 x 50µl reactions	QT615-02	€290.00
		500 x 50µl reactions	QT615-05	€550.00
Universal range of SensiMix reagents for RNA templates - suitable for probe or SYBR® Green I assays				
SensiMix One-step (For cDNA synthesis and real-time detection in a single step)				
2x SensiMix One-step, 50mM MgCl ₂ , 50x SYBR® Green I, RNase inhibitor, DEPC-treated water	2x formulation, dT based, contains internal reference	25 x 50µl reactions	QT205-01	€55.00
		250 x 50µl reactions	QT205-02	€565.00
		500 x 50µl reactions	QT205-05	€995.00
SensiMix Two-step (For cDNA synthesis and real-time detection in two separate steps)				
2x SensiMix dT, 50mM MgCl ₂ , 50x SYBR® Green I, RNase inhibitor, DEPC-treated water, Sensi reverse transcriptase, 5x reverse transcriptase buffer, 10mM dNTP mix, Oligo (dT) ₁₈ primer mix, random hexamer primer mix	2x formulation, dT based, contains internal reference	25 x 50µl reactions	QT305-01	€62.00
		250 x 50µl reactions	QT305-02	€610.00
		500 x 50µl reactions	QT305-05	€1,175.00
For High-Resolution Melt analysis (HRM™)				
SensiMix HRM				
2x SensiMix HRM, EvaGreen™ 50mM MgCl ₂	2x formulation supplied with a separate vial of EvaGreen™	250 x 25µl reactions	QT805-02	€150.00
		250 x 25µl reactions	QT805-05	€290.00

Please quote LTF/3FOR2 on your order to qualify.

Notes:

Products are for research purposes only. SensiMix™ is a trademark of Quantace Ltd. HRM™ is a trademark of Corbett Life Science. SYBR® Green I is a trademark of Molecular Probes Inc. EvaGreen™ is a trademark of Biotium Inc. Purchase of these products includes limited rights to use the supplied amount of SYBR® Green I stain patented by Molecular Probes Inc. Purchase of these products does not convey a licence to perform any patented process.

