

UV-Ersatzröhren UVP



Produkteigenschaften

- Original UVP-Ersatzröhre
- Erstausstattungsqualität
- für optimale Ergebnisse

Watt	Wellenlänge	Artikel-Nr. LTF	Artikel	Röhrenlänge* ca.
4 Watt	254 nm	110.4200	UVP 34-0003-01	135 mm
4 Watt	302 nm	110.4201	UVP 34-0071-01	135 mm
4 Watt	365 nm	110.4202	UVP 34-0005-01	135 mm
4 Watt	Weißlicht	110.4203	UVP 34-0066-01	135 mm
4 Watt	254/465 nm Multiband	110.4204	UVP 34-0004-01	135 mm
4 Watt	365 nm BLB	110.4205	UVP 34-0010-01	135 mm
6 Watt	254 nm	110.4207	UVP 34-0013-01	211 mm
6 Watt	302 nm	110.4209	UVP 34-0044-01	211 mm
6 Watt	365 nm	110.4208	UVP 34-0034-01	211 mm
6 Watt	Weißlicht	110.4210	UVP 34-0063-01	211 mm
6 Watt	254/365 nm Multiband	110.4206	UVP 34-0015-01	211 mm
8 Watt	254 nm	110.4211	UVP 34-0007-01	228 mm
8 Watt	302 nm	110.4213	UVP 34-0042-01	228 mm
8 Watt	365 nm	110.4212	UVP 34-0006-01	228 mm
8 Watt	Weißlicht	110.4214	UVP 34-0056-01	228 mm
8 Watt	365 nm BLB	110.4215	UVP 34-0031-01	228 mm
15 Watt	254 nm	110.4217	UVP 34-0008-01	435 mm
15 Watt	302 nm	110.4220	UVP 34-0039-01	435 mm
15 Watt	365 nm	110.4218	UVP 34-0009-01	435 mm
15 Watt	Weißlicht	110.4221	UVP 34-0051-01	435 mm
15 Watt	365 nm BLB	110.4219	UVP 34-0017-01	435 mm
20 Watt	254 nm	110.4222	UVP 34-0067-01	586 mm
20 Watt	365 nm BLB	110.4223	UVP 34-0068-01	586 mm
25 Watt	254 nm	110.4224	UVP 34-0073-01	435 mm
25 Watt	302 nm	110.4225	UVP 34-0072-01	435 mm
25 Watt	365 nm	110.4226	UVP 34-0060-01	435 mm
40 Watt	254 nm	110.4227	UVP 34-0067-01	1197 mm
40 Watt	365 nm	110.4229	UVP 405	1197 mm
40 Watt	365 nm BLB	110.4228	UVP 34-0068-01	1197 mm

*Röhrenlänge ohne Verbindungsstifte