



Bio-Cooler for Tubes up to -23°C

Description

Sample case for Cryo-Tubes while working or for transport. This Bio-Cooler protects your valuable enzymes and reagents from temperature fluctuations.

The Biocooler is a unique tube system consisting of an ice gel module permanently mounted in a convenient carrying case.

Includes an erasable note board on the inside of the cover for easy sample identification.

- Specially designed to protect enzymes and reagents stored down to -23°C.
- For 48 x 1.5/2.0 ml vials
- Dimensions (W x H x D): 21.5 cm x 9.5 cm x 20 cm
- Weight: 2.5 kg

