

Mastermix Tubes 5 ml

Description

5 ml tubes with screw cap.
For use with the robots PIRO, CAS-1200/4200 and other pipetting systems.

Type:

Graduated, free standing with flat bottom.
With conical shape for minimal residual.
Delivered sterile with closed cap.

Dimensions:

16 mm (tube), 20,5 mm (capl);
height 56,5 mm (open), 58 mm (closed).

Packaging Unit:

Bag with 25 pieces. Tubes with lids.



5 ml Mastermix Tubes
(Order No.: 112.3070)