

RNase FREE Spray Apex

Description

Apex RNase Free is ready-to-use and removes nuclease and nucleic acid contamination from surfaces by completely destroying the enzymes. Spray, wipe or soak labware then wipe or rinse the surface clean with RNase-free water.

Easily decontaminates a variety of surfaces and labware including your pipettor, benchtops, glassware and more. Apex RNase Free is ideal for labware which cannot be autoclaved and replace time-consuming processes such as DEPC treatments.

Order Number: 102.5000

Specification

- Easily and completely inactivates RNase and DNase contaminants on contact
- Effectively removes high levels of RNase contamination that similar products cannot
- Safe to use on pipettors, glassware, plasticware, and most metal surfaces
- Space-saving, convenient packaging
Eliminates over-spraying of decontaminants on sensitive equipment
- Large volume (475 ml)

Please Note

- Store at room temperature in a UV-resistant container.
- Skin irritant! May cause skin irritation on contact. Avoid contact with eyes and skin. Wear chemical safety goggles and chemical resistant gloves when applying.
- Eye irritant! In case of contact, rinse with plenty of water.
- Fire protection information: not flammable!

