



Mini Tankblot

Description

High-quality blotter with safety cover for the use of 4 cassettes and 8 pads.

Transfer buffer volume: 2100 ml

Suitable with vertical electrophoresis chamber H10 Mini.

Technical Data

Safety lid	Yes
Cassette capacity	4
Pad capacity	8
Electrode material	Platinum
Transfer time	60 - 120 min
Gel capacity	4 Mini gels
Blotting area (W x D)	8 x 9.4 cm
Transfer buffer volume	2100 ml
Cooling option	ja
Dimensions (W x D x H)	21 x 15 x 25 cm



Mini Tankblot



Compatible with:
Vertical gel chamber H10 Mini