



Maxi Tankblot

Description

High-quality blotter with safety cover for the use of 2 cassettes and 4 pads.

Transfer buffer volume: 6100 ml

Suitable with vertical electrophoresis chamber H20 Maxi.

Technical Data

Safety lid	Yes
Cassette capacity	2
Pad capacity	4
Electrode material	Platinum
Gel capacity	4 Maxi gels
Blotting area (W x D)	22 x 19 cm
Transfer buffer volume	6100 ml
Cooling option	ja
Dimensions (W x D x H)	32,5 x 15 x 40 cm



Maxi Tankblot



Compatible with:
Vertical gel chamber H20 Maxi