

# Gradient Mixer for Gel Electrophoresis


## Description

Device for mixing from 20 to 250 ml gradient volume.

Both mixing vessels are mounted on a stable worktop with levelling compensation.

## Applications

- Gel Electrophoresis
- Column Chromatography
- Ultracentrifugation



Gradientenmischer

Gradient Volume	Part No.
2 x 20 ml	101.0510
2 x 30 ml	101.0511
2 x 60 ml	101.0512
2 x 120 ml	101.0513
2 x 250 ml	101.0514