

# FRK Storage Systems

## Accessory Shelf for Deep Freezers

### Descriptions

The FRK storage systems are stainless steel shelves for standard containers. Due to small gaps, they optimally ensure an even and fast cooling of the complete frozen goods.

### Specifications

- Completely made of stainless steel
- For standard cryoboxes
- With handle strips
- Extra gaps for optimum, even and rapid cooling of the complete frozen product



### Technical Data

Cooling unit type	Quantity per unit and shelf type	Drawer dimensions W x D x H	Number drawers
Cooler and freezer box B 30-20	1 x B 30-RS	30 x 28 x 5,2 cm	4 drawers for a total of 16 cryoboxes
Cooler and freezer box B 35	2 x B 35-RS	30 x 28 x 5,2 cm	5 drawers for a total of 30 cryoboxes
Freezer TS	1 x TS 100-RS	30 x 28 x 5,2 cm	9 drawers for a total of 81 standard cryoboxes
Freezer TUS	1 x TS 100-RS	30 x 28 x 5,2 cm	9 drawers for a total of 81 standard cryoboxes

Delivery without cryoboxes.