



## CM-50 Centrifuge

### Description

This compact benchtop centrifuge is a unique high-speed centrifuge (15,000 RPM). The patent-pending SkySpin™ motor drive provides exceptional performance as well as quiet and smooth operation.

Part number: 109.2005

### Benefits

- Intuitive and clear interface with large color display
- secure lid lock
- Display preferences include the following:
  - Language
  - Brightness
  - Sleep mode
  - Sound volume
- Short Spin button
- High tech digital imbalance sensor with great tolerance
- Patent pending digital motor driver SkySpin™
- Advanced new chassis insure super-smooth run even when slightly imbalanced.



Extra low noise  
And heatup



Lid lock



Digital imbalance  
Sensor


### Specifications

Centrifugal speed	1.000 – 15.000
Max. centrifugal force	15.300 RCF
Timer	0.1 - 99 min
Dimensions (L x W x H)	200 x 180 x 145 mm
Power supply adapter	24 VDC, 5 A
Weight	3.1 kg


### Rotor 50.01

Universal High-Speed-Rotor for 12 test tubes 0.2 - 2 ml. Patent pending cover cap insures the optimum air flow and incredibly low sample heating.

**Max: 15.000 RPM**



Custom designed interface offers outstanding clarity and comfort.



Using the sound, brightness and language options everyone can make the machine more personal and adequate the current operation environment.