

Incubator Hood TH 15

Description

The Incubator Hood TH 15 has a temperature range from +5°C above ambient to +50°C and can be combined with all compact shakers KS 15 and Microplate Shakers TiMix 5 or used as independent unit. The height is sufficient for 1000 ml Erlenmeyer flasks.

Benefits

- Coated metal housing with 3 windows, the front door opens upwards
- Safety cut-off of the shaker if combined with a Bühler Compact Shaker
- Standard temperature controller 5°C above ambient to +50°C

Specifications

Temperature range	5°C above ambient to +50°C
Temperature accuracy	± 1% of the set value
Inner dimensions (W x D x H)	90 x 410 x 310 mm (height sufficient for 1000 ml Erlenmeyer flasks)
Temperature sensor	PT 100
Air circulation	240 m³/h
Electrical supply	230 V, 50/60 Hz
Dimensions (W x D x H)	510 x 450 x 455 mm
Weight	21 kg
Part number	104.0270
Accessories	Shelf TH 15 Additional tray for samples which have to be temperature controlled, but need no shaking. Dimensions: 378 x 318 mm.



Incubator hood TH 15



Incubator hood TH 15 with shelf