



# Rocking Shaker DRS-12

## Description

A robust, reliable and powerful rocking shaker. An attractive design, good manufacturing quality, user-friendly operation and maintenance-free nonstop operation make DRS-12 the perfect tool for countless applications in the medical, chemical and biological laboratory.

## Specifications

- Variable shaking angle from 1 to 12 degrees
- Shaking speed from 1 to 50 shakes per minute
- Custom shaking programs
- Soft run advanced motor for top performance
- Digital timer
- High-precision, operational controls
- Sound notification when operation complete
- Onboard memory for storing settings (see Page 2)



Rocking Shaker DRS-12  
Order number: 105.2006

## Technical Data

Size platform (workingspace)	257 x 307 mm
Max. load	3 kg with central arrangement, 1 kg edge
Rocking angle	1 - 12 degree
Time controller, limitless	1 - 6000 min
Rocking speed	1 - 50 rpm
Environmental temperature	+4 bis +45 °C
Relative humidity (at 20°C)	max. 80%
Power supply adapter	DC 24V, 1 A
Dimensions (L x W x H)	275 x 310 x 110 mm
Weight	3.6 kg
Order number shaker	105.2006
Order number platform	105.2120

## Display

### Programmable Mixing Modes

- Up to 4 different mixing behaviors combined in one program
- Detailed specification of each mixing step
- Various combination and cycle options



## Rocking Shaker DRS-12

### Accessories

2 platforms can be added to the rocker for maximum capacity.

Please note the maximum load (3 kg with central arrangement, 1 kg edge).



### Second platform for DRS-12

Non-slip mat inclusive  
Order number: 105.2120



### Anti-slip mat for DRS-12

Suitable for DRS-12 platforms  
Order number: 105.0524