

RollerBottle Enclosure Special Construction

Description

RollerBottle is a flexible roller mixer for bottles in various sizes with robust roller drive without up and down.

Depending on the size of the containers, individual rollers can be removed or relocated (except drive roller).

With the optional Plexiglas enclosure, the bottles or containers are additionally protected. The enclosure is specially made according to the model and your needs.

The device is then completely inside the housing and can be loaded from above. The lightweight lid is also made of Plexiglas and can be conveniently operated with a handle.

A socket for the RollerBottle is located inside on the rear panel. When the lid is opened, the power is automatically interrupted and the device stops.

Of course, the RollerBottle can still be operated independently of the housing.

A valuable enhancement to better protect samples and people present.



RollerBottle Enclosure with Roller Bottler anlaog, 650 mm



RollerBottle Enclosure