



## Mini Shaker PSU-2T

### Description

Mini Shaker PSU-2T is designed for immunoassays and provides adjustable mixing of reagents in micro-plates. The device ensures smooth movement of the platform even at low speeds.

Shaker is a compact and user-friendly device. It takes up little space on a desk and is ideal for personal use. The use of direct drive and brushless motor allows continuous mixing and trouble-free operation. Display of the device switches between time and speed readings.

Shaker can be used in cold rooms or incubators, operating at ambient temperature range +4°C to +40°C.

Item Number Minishaker PSU-2T: 105.1208



### Technical Data

Speed control range	150 – 1200 rpm
Digital time setting	1 min - 24 hrs / non-stop
Direct drive mechanism	Yes
Maximum continuous operation time	168 h
Digital setting and control of time and speed	Yes
Orbit	2 mm
Overall dimensions (W x D x H)	255 x 255 x 100 mm
Weight	2 kg
Input current/power consumption	12 V, 280 mA 3.4 W
External power supply	Input AC 100–240 V; 50/60 Hz; Output DC 12 V

### Accessories



Standard platform IPP-2



Platform IPP-4 (optional)  
Partnumber: 105.1209