



# Mini Microcentrifuge CF-5

## Description

Mini-CF-5 is a small and visually appealing microcentrifuge. It is specially designed for fast centrifugation down from the walls or lids of reaction tubes.

It offers excellent uniformity of speed control.

Perfectly suitable e.g. for the following tasks:

- Plasmid
- DNA/RNA
- Miniprep ...



Micro centrifuge CF-5 lid closed

## Specifications

- Complete set with standard rotor and strip rotor
- Easy opening by pressing the lid
- Air cooling system: prevents heat from the motor from getting inside
- Stops automatically when the lid is opened
- Quiet and low vibration
- Housing and rotor made of sturdy ABS
- Clear and transparent polycarbonate lid
- Two rotors included in delivery:  
Standard rotor for 6 x 0.5/1.5/2.0 tubes  
strip tube rotator for 2 x 0.2/0.5 ml tubes



Micro centrifuge CF-5 with standard rotor

## Technical Data

Speed setting range	max. 5.500 rpm
Max. RCF	1.775 x g
Material	Lid: PC Housing: ABS Rotor: ABS
Dimensions (W x D x H)	178 x 150 x 109 mm
Overall dimensions (W x D x H) with platform	235 x 235 x 140 mm
Weight	1.4 kg
Power	8 W
Power supply	AC 230 V 50/60 Hz
Order No.	109.0030

## Included in delivery



Strip tube rotator



Adapter for 0.5 ml tubes



Adapter for 0.2 ml tubes