








PIRO Adapters



PIRO pipetting robot provides extreme flexibility in combination with small space for various liquid handling tasks. Whether simple fluid transfer or complex PCR run, with the easy to use software your application is ready within minutes.

PIRO adapters have been created especially for PIRO pipetting robot and guarantee optimum results. Available in different well formats.

Model (well)	Description	Part number
SBS 32 well U 	32 well cooling adapter 32 x 1.5 ml / 2 ml tubes for tubes with appended caps, inclusive fixation for caps tube lifter for cooling adapter	P-S112.0326 P-S112.0054
SBS 96 well H 	96 well cooling adapter high for all kinds of non-skirted and half-skirted PCR plates, 0.2 ml PCR tubes and strip tubes	P-S112.0515
SBS 96 well L 	96 well cooling adapter flat for all kinds of non-skirted and half-skirted PCR plates, 0.2 ml PCR tubes and strip tubes	P-S112.0516
SBS 36 well 	36 well LightCycler capillary system and 0.5 ml tubes C	P-S112.0001
SBS 384 well 	for better cooling and positioning of 384 well PCR plates	P-S112.0004
SBS 96 well low universal 	96 well cooling adapter (universal) for all kind of non-skirted, half-skirted and full-skirted PCR plates, 0.2 ml PCR tubes and strips	P-S112.0063
Suitable plate racks: Adapter SBS 32 well U, SBS 96 well H and SBS 36 well L can be used with plate rack P-S112.0417 Adapter SBS 96 well L, SBS 384 well can be used with plate rack P-S112.0318 or P-S112.0319		