








PIRO Plate Racks





PIRO pipetting robot provides extreme flexibility in combination with small space for various liquid handling tasks. Whether simple fluid transfer or complex PCR run, with the easy to use software your application is ready within minutes.

PIRO plate racks have been created especially for PIRO pipetting robot and guarantee optimum results.

Type	Description	Part number
 <p>Plate Rack-1 SBS 96 well</p>	<p>For semi-skirted PCR plates and non-skirted PCR plates</p> <p>E.g.: 96 well half-skirted low profile, 96 well half-skirted, 96 well unskirted, 0,2 ml single tubes, 0,2 ml 8 strip tubes und 0,2 ml 12 strip tubes</p>	P-S112.0318
 <p>Plate Rack-1 SBS standard</p>	<p>For all kinds of full-skirted plates.</p> <p>E.g.: 384 well microtitre plate, mikrotitre plate, PCR plates</p> <p>Suited for combination with PIRO cooling adapter (P-S112.0516, P-S112.0063, P-S112.0004) and for semi-skirted PCR plates.</p>	P-S112.0319
 <p>Plate Rack-2 SBS standard</p>	<p>For all kinds of deep well plates and trays</p> <p>Suited for combination with: PIRO 96 well cooling adapter high (P-S112.0515) PIRO adapter SBS 32 well U (P-S112.0326)</p>	P-S112.0417
 <p>Plate Rack RG</p>	<p>For RotorGene GeneDisc (Corbett RotorGene and Qiagen RotorGeneQ)</p>	P-S112.0459



PIRO Plate Racks

Type	Description	Part number
Plate Rack-1 SBS SmartCycler 	For SmartCycler plates and tubes	P-S112.0466
Plate Rack-2 SBS HG 	Adapter for different third-party providers with high guides E.g.: NucleoMag	P-S112.0486
Plate Rack-1 SBS 32-well magna-pure 	For 32-well magna-pure plates of Roche	P-S112.0006
Plate Rack LC carousel V1 	For LightCycler 1.5 carousel of Roche for 32 probes	P-S112.0007
Plate Rack LC carousel V2 	For LightCycler 2.0 carousel of Roche for 32 probes	P-S112.0007