




PIRO Waste Racks



PIRO pipetting robot provides extreme flexibility in combination with small space for various liquid handling tasks. Whether simple fluid transfer or complex PCR run, with the easy to use software your application is ready within minutes.

PIRO waste racks have been created especially for PIRO pipetting robot and offer different disposal solutions.

Type	Description	Part number
Waste rack SBS 	For the collection of used tips, with special lid. Usable in combination with a waste bag	P-S112.0321
Waste box double SBS 	Waste box with extended volume (includes two SBS places)	P-S112.0005
Waste channel 	For tips with 50 µl, 200 µl and 1000 µl. Suitable for use with PIRO console and/or external waste disposal.	P-S112.0322
PIRO console 	Drawer system for tip disposal (right drawer) and storage of racks or adapters (left drawer)	P-S112.0337