




PERSONAL PRECISION

PORTABLE POWER



PORTABLE POWER SUPPLY

Get closer to your sample with a portable power supply. This compact power supply enables you to use locally available power to run your MyGo Mini.

COMPACT - CONVENIENT - COMPATIBLE

USE A BATTERY

Compatible with standard automotive batteries.

USE A CAR

Compatible with vehicle 12V socket (cigarette lighter socket).


USE A PLANE

No need to carry-on Lithium batteries.

USE IT ALL DAY

Up to 12 runs from a standard battery.

Use your power adapter to run your MyGo Mini all day.



SPECIFICATION

Input Voltage	11V to 16V DC
Input Power	100W (e.g. car battery)
Operating Conditions	0°C to 40°C, 5% to 95% RH non-condensing
Storage Conditions	-20°C to 70°C, 5% to 95% RH non-condensing
Weight	300g
Dimensions	50 x 105 x 20 mm (PSU)
Compatible with	MyGo Mini and MyGo Mini S

ORDERING DETAILS

SKU	Description
8265	MyGo Mini - PSU - 12V DC

All trademarks are the property of their respective owners. Design and specifications are subject to change without notice.

Version 001