



## IsoFreeze PCR Rack

### Description

Optimal optical control for your PCR samples!  
Made of a gel-filled middle section with temperature-dependent color changing material. So you can see at a glance whether your samples still have the right temperature.

There is a significant change in a contrasting color at 7°C, so you can use this indicator to determine when the temperature of the rack, or even a single well, rises above the normal cooling temperature. Each rack comes with a custom snap-on lid made of transparent polypropylene.



IsoFreeze PCR Rack green/yellow



IsoFreeze PCR Rack purple/pink

### Specifications

- For 96 x 0.1 or 96 x 0.2 ml tubes with snap top or screw cap.
- Maintains a sample temperature of less than 4°C for approx. 4 hours.
- Dramatic change to contrasting color at 7°C.
- Dimensions (rack + lid): 70 x 105 x 146 mm.
- 2 models available:
  - Color changes from green (cold) to yellow (warm)
  - Color changes from purpur (cold) to pink (warm)

Description	Color	Packaging unit	Part number
IsoFreeze PCR Rack Suitable for tubes: 96 x 0.1 or 96 x 0.2 ml	green/yellow	2 pieces per pack	104.0046
IsoFreeze PCR Rack Suitable for tubes: 96 x 0.1 or 96 x 0.2 ml	purple/pink	2 pieces per pack	104.0043