


## Cooling Storage Box for 15 ml / 50 ml Tubes

### Description

These polypropylene boxes are a clean and durable alternative to cardboard boxes. The clear lid ensures easy sample identification and the alpha-numeric indexing in white makes locating a tube easy.

For 15 ml and 50 ml tubes. With transparent lid. Also suitable for cold storage and temperatures down to -90°C.



### Benefits

- Robust polypropylene tube box for centrifuge tubes
- For 36 x 15 ml tubes (green box) and 16 x 50 ml tubes (blue box) with sturdy transparent lid
- Temperature range from -90°C to +121°C
- Autoclavable
- Numeric index
- 145 x 145 x 127 mm (W x D x H)

Description	Color	Packaging Unit	Part Number
Tube Storage Box for 36 x 15 ml Tubes	green	2 Boxes	113.3540
Tube Storage Box for 16 x 50 ml Tubes	blue	2 Boxes	113.3541