



## PCR CoolSafe Chambers

### Description

Solid aluminium PCR cooling units. Each tube is individually enclosed by the aluminum unit, keeping the temperature of the samples constant.




### CoolSafe for PCR Plates & 0.2 ml Tubes

Total number of holes: 106

Part No.: 104.0191

- for 0.2 ml PCR plates or tubes
- Holes:  
Center: 96 for PCR plates (8 x 12)  
Sides: 10 for 4 x 1.5 ml tubes and 6 x 0.5 ml tubes
- Hole Depth: 16 mm (5/8 Inch)
- Hole Diameter: 6.4 mm (1/4 Inch)
- Dimensions (L x B x H): 156 x 79 x 25 mm  
(6-1/8 x 3-1/8 x 1 Inches)



### CoolSafe Double for 2 PCR Plates & 0.2 ml Tubes

Total number of holes: 212

Part No.: 104.0193

- for 0.2 ml PCR plates or tubes
- Holes:  
Center: 96 for PCR plates (8 x 12)  
Sides: 10 for 4 x 1.5 ml tubes and 6 x 0.5 ml tubes
- Hole Depth: 16 mm (5/8 Inch)
- Hole Diameter: 6.4 mm (1/4 Inch)
- Dimensions (L x B x H): 156 x 158 x 25 mm  
(6-1/8 x 3-1/8 x 1 Inches)



### CoolSafe für 1.5 ml Tubes

Total number of holes: 48

Part No.: 104.0192

- for 1.5 ml tubes
- Well: 48 (4 x 12)
- Hole Depth: 35 mm (1.38 Inch)
- Hole Diameter: 6.4 mm (1/4 Inch)
- Dimensions (L x B x H): 156 x 79 x 38 mm  
(6-1/8 x 3-1/8 x 1 Inches)



### CoolSafe für 1.5 ml Tubes

Total number of holes: 80

Part No.: 104.0194

- for 1.5 ml tubes
- Well: 80 (5 x 16)
- Hole Depth: 35 mm (1.38 Inch)
- Hole Diameter: 10 mm (1/4 Inch)