

Mini UV Viewing Cabinets

Description

Mini Chromato-Vue Viewing Cabinets provide a compact darkroom environment for viewing materials. All cabinets feature a contrast control filter that absorbs UV energy to protect the eyes.

- The C-10P accommodates the rechargeable lamps which are excellent for field use applications
- Use the C-10 with compact and handheld UV lamps
- The C-10E4 is designed for use with the 4-watt EL series lamps
- The C-10E6 is designed for use with the 6-watt EL series lamps
- The C-15G cabinet combines the C-10 cabinet with the UVGL-58 multi-band shortwave/longwave handheld lamp to provide a self-contained illumination system

Applications

- Mineralogy
- Lab Fluorescence
- Sterilization
- Field Investigation

Main Features

- Each cabinet is designed to be used with specific lamp models, which fit snugly on top of the cabinet
- UV Mini viewing cabinets are constructed with a metal base and plastic top
- The bottom panel of the cabinet is removable
- A contrast control filter absorbs UV energy to protect the eyes
- Soft rubber viewport provides a comfortable viewing area with adequate space for eyeglasses

Specifications

Dimensions Interior (W x D x H): 254 x 229 x 178 mm

Dimensions Exterior (W x D x H): 267 x 229 x 305 mm

CE-certified

Models

Model	Description	Part Number
C-10	For compact/handheld lamps	110.5300
C-10P	Portable for rechargeable lamp (UVS-26P, UVSL-26P, UVL-26P)	110.5302
C-10E4	For four watt EL lamps	110.5301
C-10E6	For six watt EL lamps	110.5308
C-15G	Set: viewing cabinet C-10, 6 watt lamp UVGL-58, SW & LW 230 V	110.5303



C-10E4 Mini UV cabinet
with 4-Watt UV lamp EL-Series



C-10 Mini UV cabinet
Basic unit without lamp



C-10P Mini UV cabinet
mit UV-Handlampe



C-15G Mini UV cabinet set
with UV-Handlamp UVGL-58



Mini UV Viewing Cabinets

Models

Model Viewing Cabinet Part No.	Viewing Cabinet with Hand Lamp	Suitable UV Hand Lamp
C-10 Basic unit 110.5300		Compact and handheld lamps
C-10P (Cabinet without hand lamp) 110.5302		Portable and rechargeable lamps Models: UVS-26P, UVSL-26P, UVL-26P
C-10E4 (Cabinet without hand lamp) 110.5301		4 Watt lamps, Serie EL
C-10E6 (Cabinet without hand lamp) 110.5308		6 Watt lamps, Serie EL
C-15G (Cabinet without hand lamp) 110.5303		Set: Viewing cabinet C-10 6 Watt lamp UVGL-58, SW & LW, 230 V