



## 3UV Multi Wavelength Lamp

### Description

Innovative, versatile and patented 3UV™ lamp houses multiple wavelengths: longwave, shortwave and midrange UV in one lamp for multiple applications. Lamp is available in four, six or eight watts.

### Applications

- Inspection
- Laboratory
- Education

### Benefits

- Three wavelengths, 254/302/365 nm UV, in one lamp
- with 4, 6 or 8 watt
- Selection of each wavelength is by an easy turn of the dial located at the top of the housing
- The uniquely designed reflector surfaces mounted behind each tube provides maximum UV for fluorescence applications
- Lightweight and ergonomically designed for hand or stationary use
- Design allows the stationary use by standing the lamp on its end or side



3UV multi wavelength lamp

### Specifications

Model	Wavelength	Tubes	Dimensions (L x W x D)	Part number
3UV-34	254/302/365 nm	4 watt	241 x 114 x 76 mm	110.4000
3UV-36	254/302/365 nm	6 watt	318 x 114 x 76 mm	110.4001
3UV-38	254/302/365 nm	8 watt	394 x 114 x 76 mm	110.4002



Suitable replacement tubes are also available from us:  
info@labortechnik.com