



UV Hand Lamps

Rechargeable UV-Lamps

Description

Portable and rechargeable UV lamp. Ideal for field operations and many other mobile applications. The integrated rechargeable battery is suitable for hours of use and can be easily recharged via a power socket.

The lamps have a sturdy handle and are made of durable and lightweight plastic. Compatible with Mini Darkroom C-10P (see picture)

The model UVSL-26P can be operated simultaneously with long and short wave.

Features

- Integrated rechargeable battery, can be charged directly via the mains socket
- Equipped with two UV tubes
- 12 V adapter for car charging (optional)
- Weight: 1.8 kg
- Dimensions (L x D x H): 249 x 76 x 140 mm



Specifications

Model	UV Source	Watt	Intensity (uW/cm ²)	Part Number
UVS-26P	254 nm	2 x 6 Watt	1250	110.4059
UVSL-26P	254 / 365 nm	2 x 6 Watt	780 (SW) / 520 (LW)	110.4060
UVL-26P	365 nm BLB	2 x 6 Watt	1600	110.4058



Mini-UV-Cabinet C-10P
Part number 110.5302



12 V Adapter
Order No. 58-0127-01