



Qualitix® Reagenz Reservoirs

Description

Re-usable and autoclavable Polypropylene. Ideal for all types of applications with single and multichannel pipettes. For work in clean environment partly possible. Available with or without lid. Choose the multichannel reservoirs sterile for clean room or non-sterile.

Benefits

- Embossed graduations guarantees long lasting legibility.
- Round corners enable convenient reagent pouring.
- Bottom shape ensures bench stability.
- Hinged lid closure.
- Microplate format guarantees excellent stability.
- Stackable for space saving storage.



V-shape reservoir with lid, stackable

Model	Material	Volume	Shape	Sterilized	Packaging	Part Number
	Polypropylene	75 ml	V-shape	no	20 / box	330.01
	Polypropylene	125 ml	V-shape with lid	no	20 / box	330.01L
	Polystyrene	48 ml 8 x 6 ml	8-channel	no	30 / box	330.08
		48 ml 8 x 6 ml	8-channel	no	100 / box	330.08XL
		48 ml 12 x 4 ml	12-channel	no	30 / box	330.12
		48 ml 12 x 4 ml	12-channel	no	100 / box	330.12XL
	Polystyrene	48 ml 8 x 6 ml	8-channel	yes	10 x 1 / box	330.08.9
		48 ml 12 x 4 ml	12-channel	yes	10 x 1 / box	330.12.9