

Filter Tip Extended Length Barrier Tip BT1250

Description

Each lot undergoes stringent inspection and individual lot testing ensures the products are certified.

*Exclusive S3 - Polymer was designed to increase pipetting accuracy by virtually eliminating tip retention and sample hold-up. This prevents liquids from sticking to the tip surface but sliding residue-free over the inside surface of the tip.


All filter tips are certified free of:

- DNA
- DNase & RNase
- Endotoxine (Pyrogene)

Part number:

BT1250.N: 102.22x

BT1250 (S³)*: 102.2215



Technical Data

Graduation Marks	At 200 µl, 500 µl und 1.000 µl, s. figure
Tip Composition	Virgin Polypropylene
Tip Type Available	Natural Polypropylene (N) or S ³ -Quality*
Filter Material	High-Density Polyethylene
Pre-Sterile	Yes
Configuration	Racked (s. Page 2)
Packaging Breakdown	96 tips per rack 8 racks per pack 4 packs per case
Autoclavable	nein
Storage Conditions	Store in a clean, dry environment at room temperature 15-30°C

Pipettor Compatibility

Biohit Proline™ 1000 µl
Brand Transferpette™ 1000 µl
CLP Beta-Pette™ 1000 µl
CLP Poseidon™ 1000 µl
CLP Poseidon Electronic™ 1000 µl
Eppendorf Reference™ 1000 µl
Eppendorf Research™ 1000µl
Eppendorf Research Plus™ 1000 µl
Eppendorf Xplorer™ 1000 µl

Finnpipette™ 1000 µl & 1200 µl
Gilson Pipetman™ P1000
Hamilton™ 1000 µl
Labnet BioPette E™ 1200 µl
Nichiryo Nichipet EX™ 1000 µl
Socorex Calibri 822™ 1000 µl
VWR Ergonomic High Performance™ 1000 µl
VWR Ultra High Performance™ 1000 µl

Pipette Tips in Safety Rack

Benefits of Optimus Rack

Secured Tips

The locking clasp clicks to ensure the rack is securely closed giving you peace of mind when transporting or storing opened racks.

The rack is designed to withstand multiple drops off a lab bench, and to lock or unlock with one hand.

Stable Corners

Redesigned base has a larger footprint (SBS) allowing you to pipette of the corners without the tipping the rack.

Structural Support

A ribbing network is built into the base and snapcard to create a stable platform for multichannel pipetting, giving you secure seals on each channel.

Environmentally Conscious

The Optimus Rack is manufactured with 1/3 less plastic than leading competitor racks while providing top functionality and strength.



Sample rack

