

# Hybridization Oven HYBRI

## Description

Small footprint and big performance  
Part number: 104.0200

## Specifications

- Low Profile
- Large capacity: 6 or 12 bottles
- Plastic-coated bottle clips
- Microprocessor temperature control
- Automatic overtemperature (OTP) control
- Rapid heat up: Ambient to 65°C in 7 minutes
- Drip tray included
- Shaking platform optional
- Chamber illuminated



HYBRI with rotisserie

## Technical Data

Rotor capacity:	6 bottles 300 x 35 mm 12 bottles 150 x 35 mm
Rotation Speed:	5 RPM / 200 Volt, 6 RPM / 110 Volt
Power:	750 Watts, 220 Volt, 50/60 Hz
Temperature Range:	Ambient +5 °C to + 90 °C
Temperature Regulation:	elektronic
Temperature Display:	LED Display
Overtemperature Control:	automatic 1 °C over actual value
Dimensions Interior:	195 x 200 x 316 mm
Dimensions Exterior:	285 x 250 x 446 mm
Weight:	13 kg



Rotisserie 6/12 and drip tray