

## PIRO Pipetting Robot – Filter Tips



PIRO pipetting robot provides extreme flexibility in combination with small space for various liquid handling tasks. Whether simple fluid transfer or complex PCR run, with the easy to use software your application is ready within minutes.

PIRO filter tips have been created especially for PIRO pipetting robot and guarantee optimum results.

Volume in $\mu$ l steril / PCR-clean	Produktbild	Description	Part Number
50 $\mu$ l filter tips clear, sterilized		24 x 96 tips Total: 2304 filter tips  Certified free of: RNase, DNase, Pyrogen, DNA/RNA	P-S112.600112
200 $\mu$ l filter tips clear, sterilized		24 x 96 tips Total: 2304 filter tips  Certified free of: RNase, DNase, Pyrogen, DNA/RNA	P-S112.600113
1000 $\mu$ l filter tips clear, sterilized		24 x 96 tips Total: 2304 filter tips  Certified free of: RNase, DNase, Pyrogen, DNA/RNA	P-S112.600114
50 $\mu$ l filter tips clear, PCR-clean		24 x 96 tips Total: 2304 filter tips  Certified free of: RNase, DNase, Pyrogen, DNA/RNA	P-S112.600135
200 $\mu$ l filter tips clear, PCR-clean		24 x 96 tips Total: 2304 filter tips  Certified free of: RNase, DNase, Pyrogen, DNA/RNA	P-S112.600134
Mastermix Tubes 25 x 5 ml		One bag with 25 tubes 5 ml tubes with screw cap  Certified free of: RNase, DNase, Pyrogen, DNA/RNA	112.3070