

UV Viewing Cabinets Large / Chromato Vue

Description

The darkrooms have a removable base plate so that they can be placed over large samples or a transilluminator.

The cabinets are made of UV-stable, lightweight plastic and are CE-certified.

Models

The **C-70** Chromato-Vue Cabinet includes two each of 15-watt shortwave and longwave UV tubes mounted on one side of the cabinet.

Model **C-71** is similar to the C-70, but has four 15-watt, longwave UV tubes on one side of the unit.

The 254 nm and 365 nm UV tubes in the **C-75** are positioned on two sides of the cabinet providing uniform UV illumination for fluorescence studies.



Darkroom C-70 / C-71



Darkroom C-75

Applications

- Non-Destructive Testing
- Document Verification
- Thin Layer Chromatography

Model C-70, Order No. 110.5305

UV Source	254 nm and 365 nm
Watt	15 W (two tubes each wavelength)
Light Sources	Overhead white light
Dimensions Interior (W x D x H)	432 x 356 x 127 mm
Dimensions Exterior (W x D x H)	594 x 400 x 330 mm
Weight	10.3 kg



WARNING:

Do not expose eyes and skin to short wave ultraviolet radiation without protection, as the radiation may be harmful to unprotected eyes and skin.

Model C-71 Order No. 110.5306

UV Source	365 nm
Watt	15 W (four tubes)
Light Sources	Overhead white light
Dimensions Interior (W x D x H)	432 x 356 x 127 mm
Dimensions Exterior (W x D x H)	594 x 400 x 330 mm
Weight	10.3 kg

Model C-75 Order No. 110.5307

UV Source	254 nm und 365 nm
Watt	15 W (two tubes each wavelength)
Dimensions Interior (W x D x H)	432 x 356 x 167 mm
Dimensions Exterior (W x D x H)	602 x 419 x 400 mm
Weight	8 kg