



# UV Viewing Cabinet C-65 / Chromat Vue

## Description

Multifunctional UVP Mid Size Chromato-Vue Viewing Cabinet provides a darkroom environment for viewing materials.

The C-65 Chromato-Vue Cabinet is designed to accommodate one or two EL Series 8-watt lamps (ordered separately, see page 2)



C-65 UV Cabinet

## Benefits

- Use of two lamps illuminates the cabinet's interior from two sides of the cabinet
- Lamps are removable from the cabinet for use in other applications
- The bottom panel of the cabinet is removable
- Soft rubber viewport and contrast control filter that absorbs UV energy to protect the eyes
- Lightweight access curtain blocks external light

## Applications

- Thin Layer Chromatography
- Non Destructive Testing
- Inspection
- Mineralogy
- Fluorescence



UVL-28LG Hand lamp suitable for UV viewing cabinet C-65

## Specifications

C-65 Light sources	White Light UV Lamps - EL series lamps ordered separately
Power switches	Located on top of the unit
Viewport	With UV blocking filter
Dimensions Interior (W x D x H)	432 x 356 x 127mm
Dimensions Exterior (W x D x H)	502 x 419 x 337 mm
Netzspannung	230 V
Certification	CE

Type	Description	Part No.
C-65	Cabinet, Base with White Light (no UV source)	110.5304
UVL Hand Lamp	Located on top of the unit	see page 2

## Handlamps suitable with UV Viewing Cabinet C-65

Lamps come with a universal plug which connects directly to the C-65 Cabinet, plus a supplementary power cord to plug into a line current for operating the lamp in other applications.



Model	Wavelength	Intensity (uW/cm <sup>2</sup> )	Part No.
UVL-18	8 Watt, 365 nm / Weißlicht	500	110.4016
UVL-28	8 Watt, 365 nm	2300	110.4019
UVLM-28	8 Watt, 2 UV, 365 nm / 302 nm	500 (LW) / 640 (MW)	110.4022
UVLMS-38	8 Watt, 3 UV, 365 / 302 / 254 nm	500 (LW) / 640 (MW) / 630 (SW)	110.4024
UVLS-28	8 Watt, 2 UV, 365 nm / 254 nm	500 (LW) / 630 (SW)	110.4023
UVM-18	8 Watt, 302 nm / Weißlicht	640	110.4017
UVM-28	8 Watt, 302 nm	1300	110.4020
UVS-18	8 Watt, 254 nm / Weißlicht	630	110.4018
UVS-28	8 Watt, 254 nm	1300	110.4021