

MyGo Quant – Performance Verification Kit

Description

The MyGo Quant is a performance verification kit for use on MyGo instruments. It is easy to use and requires no extra reagents.

Quant-Kit verifies the reliability of qPCR applications on your MyGo cycler.

Shipping and Storage

Unopened and stored at room temperature, the MyGo Quant-Kit has a shelf life of 12 months. If the kit has been opened, it must be stored at 4°C and used within 12 hours.




Content

Component	Cap Colour
1 x MyGo Buffer	Clear
1 x MyGo Probes Master	Black (glass vial)
1 x DNA Dilution 1	Red
1 x DNA Dilution 2	Yellow
1 x DNA Dilution 3	Green
1 x DNA Dilution 4	Blue
1 x No Template Control (NTC)	Black

Product	Order No.
MyGo Quant performance verification kit (48 reactions)	102.1515*
MyGo Quant performance verification kit (48 reactions) pack of 10	102.1520

*MyGo Quant kit provides a total of 48 reactions, sufficient to check the performance of one MyGo PRO/ESR and one MyGo Mini.



The results shown are typical for MyGo cyclers. If you are unsure about your results or need support, please contact LTF Labortechnik. Please note the batch number of the test kit in your experiment under notes for our support.