



# Stainless Steel Freezer-Rack

## Description

Horizontal stainless steel racks with fully extendable drawers.

Suitable for 15 ml and 50 ml tubes.

All tubes are stored upright and fit all 15 ml and 50 ml tubes.

The drawer can be pulled out completely for free access to all samples.

All tubes are organized in a grid system that ensures even circulation.



## Models

Tube volume	15 ml	15 ml	50 ml	50 ml
Shelf size	419 x 140 x 252 mm	559 x 140 x 252 mm	559 x 140 x 252 mm	419 x 140 x 252 mm
Capacity total	120 Tubes	160 Tubes	78 Tubes	60 Tubes
Part no.	113.3515	113.3520	113.3530	113.3620