

# Cold Box B 30-20

## Description

The cooling and freezing box B 30-20 is suitable as a tabletop unit directly for installation at the workplace. The units are stackable. Up to two devices can be placed on top of each other.

The model 30 with a useful capacity of 30 liters has an insulating glass pane in the door, through which the through which the refrigerated goods can be observed.

## Specifications

- Housing made entirely of stainless steel
- Cooling unit: low-noise, fully hermetic, air-cooled, low-maintenance
- Electronic temperature control with LED display
- Interior without fixtures allows easy cleaning and disinfection
- Cable gland ( $\varnothing$  10 mm), e.g. for temperature monitoring with external sensors
- Lockable door lever
- Door with insulating glass
  - Double magnetic seal with additional lip
  - Door frame heating
- Stackable units (max. 3 units)
- Impact protection rings all round at top and bottom



Cold Box B 30-20

## Option

Recirculation fan in the interior for better temperature distribution.

## Accessories

### Stainless Steel Shelf B 30-RS

- For 16 standard cryoboxes (135 x 135 x 52 mm)
- Four drawers for four standard cryoboxes each
- Footprint per drawer (W x D x H):  
300 x 280 x 52 mm
- Extra spaces for optimal, uniform and fast cooling of the complete of the complete goods



Stainless Steel Shelf B 30-RS

## Cold Box B 30-20



### Technical Data

Temperature range	-20°C to +10°C
Control accuracy	±1 K
Capacity	30 litres
Interior dimensions (W x D x H)	360 x 350 x 230 mm
Exterior dimensions (W x D x H)	530 x 700 x 460 mm
Weight	42 kg
Ambient temperature (min/max)	+12°C to +30°C
Electrical connection	230 V / 50 Hz
Order No. Cold Box B 30-20	104.5100
Order No. Shelf	104.5110