

Pen-Ray® Field-Pack

Description

Field-Pack combines a Pen-Ray® UV light source and a power supply unit (9V battery) in a compact, portable housing. The small Pen-Ray® lamps are very lightweight and are perfect for taking with you on field studies.

Applications

- Field use inspections
- Spectroscopy
- Fluorometry
- Wavelength calibration

Specifications


- Field Packs weigh only about 140 and fit in a pocket or can be worn on a belt with the clip attached to the case.
- Field Packs are available with a choice of customisation options or in three standard configurations.
- Model 95-0256-03 contains two interchangeable light sources (365 nm and 254 nm). The lamp connector is located on the outside of the power supply (this configuration only).
- Other lamp configurations are equipped with one light source and permanently connected to the power supply unit.



Pen-Ray® Field-Pack

Technical Data

Order No.	Lighted Length	Ø Quartz	Overall Length	Ø Handle	Wavelength
95-0256-01	44.5 mm	6.5 mm	108 mm	9.5 mm	302 nm
95-0256-02	44.5 mm	6.5 mm	108 mm	9.5 mm	254 nm
95-0256-03	44.5 mm	6.5 mm	108 mm	9.5 mm	365 nm and 254 nm



Construction of Pen-Ray® lamps