

RCS Incubator Hood

Description

This portable RCS incubator hood is ideal for all small-scale incubator applications. Suitable for:

- Shaker
- Magnetic Stirrer ...
- Temperature sensitive equipments with maximum dimensions of 69 x 65,5 x 61 cm (W x D x H)

Specifications

- Controlled by microprocessors
- Over temperature protection
- Digital display of actual and target values
- Temperature accuracy of $\pm 0,2^{\circ}\text{C}$

Benefits

- A big front panel with self adjustment makes easy to handle and control the installed equipment.
- Acrylic glass is the best material to protect the temperature.
- Constant temperature from ambient up to $+65^{\circ}\text{C}$
- Accuracy $\pm 0,2^{\circ}\text{C}$.
- As an option, we offer a cooling coil system to drop the chamber temperature down to $+8^{\circ}\text{C}$ in reference to the ambient temperature of 20°C (Model UH 4/14).



Option

Cooling Function with Recirculating Chiller Model UH4/14 D

Cooling for Incubators
Cooling capacity up to 250 W
Item No.: 104.0050



UH4/14 Recirculating Chiller

Technical Data

Model	RCS-550	RCS-650	RCS-750
Internal Dimensions (W x D x H)	49 x 65.5 x 61 cm	59 x 65.5 x 61 cm	69 x 65.5 x 61 cm
Width of Front Space	56 + 10.5 cm	66 + 10.5 cm	76 + 10.5 cm
Airflow Capacity	260 m ³ /h	260 m ³ /h	260 m ³ /h
Power Source	230 V, 50/60 Hz	230 V, 50/60 Hz	230 V, 50/60 Hz
Item No.	104.0250	104.0251	104.0252