

MULTIDIZER HM-4000

Description

MULTIDIZER Oven is designed with two separate incubation chambers that work as separate, independent units and can be used for a variety of applications.

The two chambers allow for simultaneous hybridization and blotting procedures requiring different sized bottles, motions, or temperature settings.

The upper chamber is equipped with a shaker tray that can be exchanged with an acrylic carousel or an orbital tray that are available as optional accessories. The removable protective trays allow for easy cleanup of spilled media. Through the use of accessories, the UVP Multidizer Oven offers several motion functions: shaking, rolling, orbital rocking and rotating all in one unit.

The unit features a precise microprocessor for control of rotation speeds and temperature control. An internal circulating fan maintains temperature uniformity, and a large LED display shows the current chamber temperature. Internal stainless steel construction enables the environment to be easily kept clean.

MULTIDIZER oven accommodates multiple bottle sizes, arrangements, and offset positions to meet a variety of mixing requirements. The rotary wheel holds twenty 35 x 150 mm bottles or 35 x 100 mm bottles, ten 35 x 300 mm bottles, or a combination of the various sizes.

Benefits

- LED display for easy operation
- Easy to clean interior made of stainless steel
- Two independent incubation chambers
- Glass door for easy viewing of samples and blocking of β -emission from the hybridization chamber
- Microprocessor controlled incubation temperature
- Accurate temperature and rotation speed control
- Carousel for multiple bottle sizes and off-set bottle positioning



MULTIDIZER HM-4000

Optional Accessories



Orbital shaker
30 U/min, angle 2°, max. load: 1.36 kg



Shaker with with back and forth movement
254 x 152 mm

MULTIDIZER HM-4000



Technical Data

Order No. MULTIDIZER	104.1603	
Total number of chambers	2	
Chambers	Lower chamber	Upper chamber
Min. temperature	10 °C above RT	
Max. temperature	99,9 °C	80 °C
Temperature accuracy	± 0.3° C to 68° C	
Temperature	± 0.1° C to 68° C	
Rotor speed	10 – 18 rpm	12 – 20 rpm
Bottles (30 cm)	10	4
Bottles (10 cm)	20	8
Bottles (15 cm)	20	8
Shaker tray	-	54 – 106 rpm
Dimensions (W x D x H)	450 x 380 x 720 mm	
Dimensions Hybridization chambers	356 x 273 x 273 mm	356 x 273 x 273 mm
Weight	34 kg	
Electrical	230 V, 50/60 Hz	
Heating element	1.250 W	500 W