

96 Well - SBS Aluminium Rack

Description

Sturdy aluminium rack in SBS format.
For 96 x 0.2 ml or 0.5 ml PCR tubes and
96-well PCR plates.
Perfect for keeping samples cool for longer during
your PCR setup.

Modern and high-quality design in shiny metallic
blue. Suitable for all SBS standard formats.

Four small rubber pads on the underside ensure
noiseless handling on smooth surfaces. Gentle on
sensitive surfaces.

The following tubes are compatible:

- 0.2 ml and 0.5 ml single tubes
- 0.2 ml strip tubes
- 96-well PCR plates

Order No.: 112.3023PCR



96 well aluminium rack, SBS format