

IsoFreeze PCR-Rack, SBS

Description

Perfect optical control for cooled samples!
Special SBS format for automated systems.
The middle part of the rack is filled with a gel.
This gel is able to change colour depending on temperature.
So you can immediately see whether your samples still have the right temperature.

A significant change in color at 7°C indicates that the temperature of the rack or even a single well is rising above the normal cooling temperature.



IsoFreeze SBS rack green/yellow

Specifications

- SBS footprint, without frames
- Ideal for automated systems, without lid
- Two models:
Color change from purple (cold) to pink (warm)
- Color change from green (cold) to yellow (warm)
- Keeps a sample temperature of less than 4°C for up to 4 hours
- Significant color change at 7°C
- 96 x 0,1 and 96 x 0,2 ml tubes
- For single tubes, strip tubes or plates
- Dimensions (W x D x H): 128 x 85 x 26 mm



IsoFreeze SBS rack purple/pink

Description	Color	Packaging Unit	Item No.
IsoFreeze SBS rack Suitable tube sizes: 96 x 0,1 oder 96 x 0,2 ml	green/yellow	2 pieces / pack	104.0048 (27-401S)
IsoFreeze SBS rack Suitable tube sizes: 96 x 0,1 oder 96 x 0,2 ml	purple/pink	2 pieces / pack	104.0045 (27-400S)