

Polypropylene Transport and Storage Box 100 Places

Description

Polypropylene storage boxes are a clean and durable alternative to cardboard boxes.

Storage box made of polypropylene with a fixed subdivision of 10 x 10 spaces for storing up to 100 tubes with a maximum diameter of 13 mm.

Fits into existing standard storage modules and offers space for all 1.5 ml to 2.0 ml microcentrifuge tubes and standard cryovials. Suitable for use in refrigerators and freezers down to a maximum of -90°C.



Specifications

- 100 places with transparent lid
- 9 x 9 format
- Lockable hinged lid for stable stacking
- Bright, attractive colors
- Stores 1.5 or 2.0 ml microcentrifuge tubes or screw cap tubes (max. Ø 13 mm)
- Temperature range: -90°C to -121°C
- Autoclavable
- Embossed reference grid in the base
- Dimensions:
151 x 151 x 54 mm (L x W x H)
- Can be used with a rack for freezers specially developed for this storage box

Accessories

Stainless steel racks for chest freezers. Specially developed for use with the 100-well storage boxes. They have a handle at the top and a latch that prevents the boxes from sliding out when the rack is inserted into or removed from the freezer.



Stainless Steel rack vertical for 11 boxes

Polypropylene Storage Box 100 Places






Figure	Colour	Order No.	Description Packaging unit
	Neon assorted	33-206AX	Polypropylene storage boxes with 100 places Suitable for 1.5 – 2 ml tubes same-colored lid, attached Packaging unit: 5 boxes
	natural	33-206X	Polypropylene storage boxes with 100 places Suitable for 1.5 – 2 ml tubes same-colored lid, attached Packaging unit: 5 boxes natural
	blue	33-206BX	Polypropylene storage boxes with 100 places Suitable for 1.5 – 2 ml tubes same-colored lid, attached Packaging unit: 5 boxes blue
	green	33-206GX	Polypropylene storage boxes with 100 places Suitable for 1.5 – 2 ml tubes same-colored lid, attached Packaging unit: 5 boxes green
	orange	33-206OX	Polypropylene storage boxes with 100 places Suitable for 1.5 – 2 ml tubes same-colored lid, attached Packaging unit: 5 boxes orange
	pink	33-206PX	Polypropylene storage boxes with 100 places Suitable for 1.5 – 2 ml tubes same-colored lid, attached Packaging unit: 5 boxes pink
	yellow	33-206YX	Polypropylene storage boxes with 100 places Suitable for 1.5 – 2 ml tubes same-colored lid, attached Packaging unit: 5 boxes yellow