

Polypropylene Transport and Storage Box 50 Places

Description

Polypropylene storage boxes are a clean and durable alternative to cardboard boxes.

Storage box made of polypropylene with a fixed partition of 5 x 10 spaces for storing up to 50 containers with a maximum diameter of 13 mm.




Specifications

- Storage box made of polypropylene with fixed
- 5 x 10 subdivision for the storage of up to 50 containers with a maximum diameter of 13 mm
- For 1.5 / 2.0 ml microcentrifuge tubes or screw cap tubes
- Lockable hinged lid with positioning lugs for stable stacking
- Embossed reference grid in the lower part.
- For the temperature range -90 °C to 121 °C
- Autoclavable
- Dimensions (W x D x H): 140 x 80 x 55 mm

Polypropylene Storage Box 100 Places



Fig. color example*	Color	Description Packaging unit
	Neon assorted	Polypropylene Storage Boxes with 50 places Suitable for 1.5 – 2.0 ml tubes Same-colored lid attached VE: 5 boxes
	natural	Polypropylene Storage Boxes with 50 places Suitable for 1.5 – 2.0 ml tubes Same-colored lid attached VE: 5 boxes natural
	blue	Polypropylene Storage Boxes with 50 places Suitable for 1.5 – 2.0 ml tubes Same-colored lid attached VE: 5 boxes blue
	green	Polypropylene Storage Boxes with 50 places Suitable for 1.5 – 2.0 ml tubes Same-colored lid attached VE: 5 boxes green
	orange	Polypropylene Storage Boxes with 50 places Suitable for 1.5 – 2.0 ml tubes Same-colored lid attached VE: 5 boxes orange
	pink	Polypropylene Storage Boxes with 50 places Suitable for 1.5 – 2.0 ml tubes Same-colored lid attached VE: 5 boxes pink
	yellow	Polypropylene Storage Boxes with 50 places Suitable for 1.5 – 2.0 ml tubes Same-colored lid attached VE: 5 boxes yellow

*Only color green original image, other figures only color example