

Polypropylene PCR Rack 96 Well

Description

Polypropylene rack for sample storage or for use as a work rack. For a variety of applications for sample processing before and after PCR.

The standard formatted rack is supplied with a transparent snap-on lid.

They are guaranteed to be dimensionally stable and can withstand temperatures down to -90°C.

Specifications

- For all standard PCR strip tubes and plates (96 x 0.2 ml)
- Lid included
- Stackable with and without lid
- Durable polypropylene for temperatures down to -90°C
- Dimensions with lid: 138 x 98 x 39 mm (L x W x H)
- Attention: not autoclavable!



96 well PCR rack

Figure	Color	Description Packaging Unit
	5 colors assorted Colors as shown	96-Well PCR-Rack Inklusive Deckel VE: 5 Boxen