

Vertical Stainless Steel Shelf for Standard Tube Boxes

Description

Stainless steel shelves for storing sample boxes in freezers.

A practical handle on the top of the shelf makes it easy to pull out the samples. An additional latch prevents the boxes from slipping out during removal.

Ideal for storing sample boxes in freezers at -20°C and -90°C.

Suitable for use with standard storage boxes with maximum dimensions of 136 x 136 x 55 mm.



Vertical stainless steel shelf,
Example figure with 11 compartments.

Models

Capacity	Dimensions
8 compartments	147 × 159 × 479 mm
9 compartments	147 × 159 × 538 mm
10 compartments	147 × 159 × 596 mm
11 compartments	147 × 159 × 655 mm
12 compartments	147 × 159 × 714 mm