

UV Handheld Lamps with Grip

Description

These UV hand lamps are characterized by their ergonomic design and are available in short, medium and long wave.

A special feature is the model UVGL-58 with two tubes in different wavelengths. These can be selected separately with a simple switch.



UVGL-55 Handl lamp with grip



Specifications

- Light and stable plastic housing
- Ergonomic handle
- 6 Watt tubes
- Weight 0.9 kg
- Combined short-/longwave version: selection possible by means of a simple click lock
- Lamp size (L x W x H): 378 x 81 x 64 mm
- Optional: matching lamp stand J-129
Size (H x L x W): (273 x 330 x 76 mm)

Models

Model	Wavelength	Tubes	Description	Item No.
UVG-54	254 nm	6 Watt	UV lamp with grip	110.4038
UVGL-55	254 / 365 nm	6 Watt	UV Lamp with grip 1 Multi-band split tube	110.4060
UVGL-58	254 / 365 nm	6 Watt	UV lamp with grip, two tubes with different wavelengths	110.4034
UVL-56	365 nm	6 Watt	UV lamp with grip	110.4036
UVM-57	302 nm	6 Watt	UV lamp with grip	110.4035



Caution!

Exposure to UV light is harmful and can cause damage to unprotected eyes and skin. The operator should wear a closed lab coat, UV certified safety glasses and gloves, which should overlap the lab coat cuffs.

Zubehör



UV Viewing Cabinet C-10
Order No.: 110.5300



Lamp stand J-129
Order No.: 110.4033