

96 x 0.1 ml Platte Low Profile

Beschreibung

Diese Platten sind speziell für PCR- und qPCR-Anwendungen konzipiert die Low Profile (LP) Produkte akzeptieren. Sie sind fest, flach und vollständig mit einem Rand versehen (full skirted). Die Platten können in einer Reihe von 0,1 ml und 0,2-ml Blöcken, PCR- und Real-Time-Thermocyclern verwendet werden.

Die Bedeutung der Einwegartikel für eine (q)PCR-Reaktion wird oft unterschätzt. Diese PCR-Platten sind hochwertige (q)PCR-Gefäße und tragen erheblich zu reproduzierbaren und verlässlichen Ergebnissen bei.

Profitieren Sie von den Vorteilen und der Benutzerfreundlichkeit von (q)PCR-Platten höchster Qualität. Sie weisen immer die gleiche Zusammensetzung des Rohmaterials und eine gleichmäßige Wandstärke auf. Variablen bei der Einheitlichkeit und der Materialzusammensetzung können damit ausgeschlossen werden.

Spezifikationen

- Robust, stapelbar und roboter-freundlich
- kein Verschachteln beim Stapeln
- Farbcodierung mithilfe eines Pre-Post Raster möglich
- Verschluss mit Strip-Caps oder Opti-Seals™ (Bestellnummer Opti-Seals: 102.0080)
- Garantierte und getestete Abwesenheit von nachweisbarer DNase/RNase, DNA/RNA, Pyrogenen und Metallen.
- ohne Autoklavieren verwendbar

Kompatible PCR-Cycler-Modelle finden Sie auf den Seiten 2-3.



96 x 0,1 ml Platte, weiß

Technische Daten

Bestellnummer	102.0079
Verpackungseinheit	25 Platten/Packung
Format	96 x 0,1 ml (Thin wall) Full skirted, für Robotic geeignet
Low Profile	Höhe: 15,6 – 16,6 Für die meisten PCR Cycler geeignet; mit verbesserten Verdunstungseigenschaften
Farbe	weiß: für die höchsten Signal-Rausch-Verhältnisse
qPCR Grad 2	Überlegenes Signal-Rausch-Verhältnis durch weiße Farbe
Verdunstungs-kategorie 1	Der Verlust des PCR-Volumens beträgt weniger als 3% PCR ist bereits mit 5µl möglich
Ausrüstung	Sehr geringe Haftung Speziell behandelt, um die Proteinbindung zu reduzieren.
Garantiert frei von	DNA/RNA, DNase/RNase, Pyrogenen und Metallen
Steril	nein
Autoklavierbar	ja

96 x 0.1 ml Platte, Low-Profile



Informationen zur Geräteunterstützung

Analytik-Jena	Alpha SC	96 x 0.2ml
Analytik-Jena	FlexCycler	Monoblock 96 x 0.2ml
Analytik-Jena - Biometra	TOne 96	96 x 0.2ml
Analytik-Jena - Biometra	TOne 96G	96 x 0.2ml
Axygen®	MaxyGene™ II	96 x 0.2ml
Bioer	GenePro B-96G	96 x 0.2ml
Bioer	LifePro TC-96/G/H(b)A	96 x 0.2ml
Bioer	Line-Gene 9600	96 x 0.2ml
Bio-Rad Laboratories	C1000	96-Deep Well
Bio-Rad Laboratories	C1000	96-Well Fast
Bio-Rad Laboratories	C1000 Touch	96-Deep Well
Bio-Rad Laboratories	C1000 Touch	96-Well Fast
Bio-Rad Laboratories	CFX96 Connect Real-Time PCR Detection System	96 x 0.2ml
Bio-Rad Laboratories	CFX96 Connect Real-Time PCR Detection System	96-Deep Well
Bio-Rad Laboratories	CFX96 Opus	96 x 0.2ml
Bio-Rad Laboratories	CFX96 Real-Time PCR Detection System	96 x 0.2ml
Bio-Rad Laboratories	CFX96 Real-Time PCR Detection System	96-Deep Well
Bio-Rad Laboratories	DNA Engine	96v Alpha unit 96 x 0.2ml
Bio-Rad Laboratories	DNA Engine	Chromo4 96 x 0.2ml
Bio-Rad Laboratories	DNA Engine Dyad	96v Alpha unit 96 x 0.2ml
Bio-Rad Laboratories	DNA Engine Dyad	Chromo4 96 x 0.2ml
Bio-Rad Laboratories	DNA Engine Dyad Disciple	96v Alpha unit 96 x 0.2ml
Bio-Rad Laboratories	DNA Engine Dyad Disciple	Chromo4 96 x 0.2ml
Bio-Rad Laboratories	DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml
Bio-Rad Laboratories	DNA Engine Tetrad 2	Chromo4 96 x 0.2ml
Bio-Rad Laboratories	PTC-100	96 x 0.2ml
Bio-Rad Laboratories	S1000	96-Deep Well
Bio-Rad Laboratories	S1000	96-Well Fast
Blue-Ray Biotech	TurboCycler TCST-9610	96 x 0.2ml
Blue-Ray Biotech	TurboCycler TCST-9620	96 x 0.2ml
Bulldog Bio	GenePro™ Thermal Cycler	96 x 0.2ml
Bulldog Bio	GeneTouch™ Thermal Cycler	96 x 0.2ml
Bulldog Bio	LifePro™ Thermal Cycler	96 x 0.2ml
Bulldog Bio	LifeTouch™ Thermal Cycler	96 x 0.2ml
Corbett Research	Palm-Cycler	96 x 0.2ml
eENZYME	GenePro Thermal Cycler	96 x 0.2ml
eENZYME	LifeExpress B-96D	96 x 0.2ml
eENZYME	LifeExpress TC-96	96 x 0.2ml
Eppendorf	Mastercycler	96 x 0.2ml
Eppendorf	Mastercycler	96-77 96 x 0.2ml
Eppendorf	Mastercycler ep gradient	96 x 0.2ml
Eppendorf	Mastercycler ep gradient Silver	96 x 0.2ml
Eppendorf	Mastercycler pro	96 x 0.2ml
Eppendorf	Mastercycler pro S	96 x 0.2ml
Eppendorf	Mastercycler X50a	96 x 0.2ml
Eppendorf	Mastercycler X50i	96 x 0.2ml
Eppendorf	Mastercycler X50l	96 x 0.2ml
Eppendorf	Mastercycler X50s	96 x 0.2ml
Eppendorf	Mastercycler X50t	384 Microplate
Eppendorf	Thermomixer C	96 x 0.2ml
Esco	Swift™ Maxi	96 x 0.2ml
Esco	Swift™ Maxi	96 x 0.2ml gradient
Esco	Swift™ MaxPro	96 x 0.2ml

96 x 0.1 ml Platte, Low-Profile



Informationen zur Geräteunterstützung

GE Healthcare	MegaBACE 1000	96 x 0.1ml
GE Healthcare	MegaBACE 4000	96 x 0.1ml
GE Healthcare	MegaBACE 500	96 x 0.1ml
Hain-LIFESCIENCE	GTQ-Cycler 96	96 x 0.2ml
Kyratec	SuperCycler	96 x 0.2ml
LongGene Scientific Instruments Co.,Ltd.	A100	96 x 0.2ml
LongGene Scientific Instruments Co.,Ltd.	A200	96-77 77 x 0.5ml
LongGene Scientific Instruments Co.,Ltd.	A200	96-77 96 x 0.2ml
MJ Research Inc.	PTC-100	96 x 0.2ml
MJ Research Inc.	PTC-100	96 x 0.2ml AgV
MJ Research Inc.	PTC-200 DNA Engine	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-200 DNA Engine	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-200 DNA Engine	Chromo4 96 x 0.2ml
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-220 DNA Engine Dyad	Chromo4 96 x 0.2ml
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-221 DNA Engine Disciple	Chromo4 96 x 0.2ml
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-225 DNA Engine Tetrad	Chromo4 96 x 0.2ml
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Moto Alpha Unit 96 x 0.2ml
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	Chromo4 96 x 0.2ml
PCRMax	AlphaCycler AC-1	96 x 0.2ml
PCRMax	AlphaCycler AC-2	96 x 0.2ml
PCRMax	AlphaCycler AC-4	96 x 0.2ml
Qiagen	QIAamplifier	96 x 0.2ml
Qiagen	QIAquant 96 2plex	96 x 0.2ml
Qiagen	QIAquant 96 5plex	96 x 0.2ml
Techne	TC-4000	96 x 0.2ml
Techne	TC-5000	96 x 0.2ml
Techne	TC-Plus	96 x 0.2ml
Thermo Hybaid	MultiBlock System	96 x 0.2ml
ThermoFisher Scientific	Arktik Thermal Cycler	96 x 0.2ml
VWR	DuoCycler	96 x 0.2ml
VWR	DuoCycler fast	96 x 0.2ml
VWR	QuattroCycler	96 x 0.2ml
VWR	QuattroCycler fast	96 x 0.2ml
VWR	UnoCycler	96 x 0.2ml
VWR	UnoCycler fast	96 x 0.2ml
Whatman Biometra	T1 Thermocycler	96 x 0.2ml
Whatman Biometra	TGradient	96 x 0.2ml
Whatman Biometra	TProfessional	96 x 0.2ml
Whatman Biometra	TProfessional Standard Gradient	96 x 0.2ml
Whatman Biometra	TRobot	96 x 0.2ml