

UVC Disinfection Lamps AR Series 300–1500

Description

AR units are used to sterilise room air and surfaces and are available in various lengths. They were designed for use without the presence of persons.


Air, surfaces and work areas are sterilized to eliminate germs, viruses, yeasts and mold spores using UV-C light.

Applications

- Laboratories
- Sterile workbenches
- Incubators
- Safety workbenches
- Work surfaces
- Production surfaces
- Storage rooms

Specifications

Series AR 300–1500 UVC Disinfection Lamp



Technical Data

Model	AR300	AR400	AR600	AR810	AR1000	AR1200	AR1500
Part No.	110.4999	110.5500	110.5501	110.5502	110.5503	110.5504	110.5505
Dimensions mm (LxWxH)	390x107x52	490x107x					