

Reagent Reservoirs Biotix

Specifications

- Biotix reservoirs are manufactured under stringent controls in an ISO 9001 certified facility.
- Biotix reagent reservoirs are certified free of DNase, RNase, and endotoxins (pyrogens).
- Every lot is tested for engineering compliance and accuracy.
- Designed for use in multichannel pipetting applications and fully compatible with 12-channel pipettes. Also available in an SBS footprint for automation platforms.
- Angled sides and a V-shaped bottom allow for maximum reagent recovery.
- Manufactured with internal volume graduation lines and pour-spouts on the corners.
- Available in polypropylene (PET) or polystyrene (PS).



50 ml reagent reservoir, polypropylene (PET)

Technical Data

Model	Volume	Description	Pre-Sterile	Quantity	Order No.
	25 ml	Material: Polystyrene (PS) Multi-channel 12-channel pipettes Angeled sides V-shaped bottom Internal volume graduation lines por-spouts on the corners SBS compliant	No	Bulk 100 reservoirs total	SR-0025-BWM
			Yes	1 reservoir/pack 25 in total	SR-0025-1SWM
				5 reservoirs/pack 100 in total	SR-0025-5SWM
	25 ml 3 Well	Material: Polystyrene (PS) Multi-channel 12-channel pipettes Angeled sides V-shaped bottom Internal volume graduation lines por-spouts on the corners SBS compliant	No	Bulk 100 reservoirs in total	SR-0025-B3WM
			Yes	1 reservoir/pack 25 in total	SR-0025-13SWM
				5 reservoirs/pack 100 in total	SR-0025-53SWM

VBM_016_E_03/2026 Subject to technical changes and errors.

Reagenz-Reservoire Biotix

Model	Volume	Description	Pre-Sterile	Quantity	Order No.
	50 ml	Material: Polypropylene (PET) Multi-channel 12-channel pipettes Angeled sides V-shaped bottom Internal volume graduation lines por-spouts on the corners SBS compliant	No	Bulk 100 reservoirs in total	SR-0050-BC
			Yes	1 reservoir/pack 100 in total	SR-0050-1SC
				5 reservoirs/pack 200 in total	SR-0050-5SC
	55 ml	Material: Polystyrene (PS) Multi-channel 12-channel pipettes Angeled sides V-shaped bottom Internal volume graduation lines por-spouts on the corners SBS compliant	No	Bulk 100 reservoirs in total	SR-0055-BNM
			Yes	1 Rrservoir/pack 25 in total	SR-0055-1SNM
				5 reservoirs/pack 100 in total	SR-0055-5SNM
	100 ml	Material: Polypropylene (PET) Multi-channel 12-channel pipettes Angeled sides V-shaped bottom Internal volume graduation lines por-spouts on the corners SBS compliant	No	Bulk 100 reservoirs in total	SR-0100-BC
			Yes	1 reservoir/pack 100 in total	SR-0100-1SC
				5 reservoirs/unit 200 in total	SR-0100-5SC
	100 ml	Material: Polystyrene (PS) Multi-channel 12-channel pipettes Angeled sides V-shaped bottom Internal volume graduation lines por-spouts on the corners SBS compliant	No	Bulk 100 reservoirs in total	SR-0100-BNM
			Yes	1 reservoir/pack 100 in total	SR-0100-5SNM
				5 reservoirs/pack 200 in total	SR-0100-1SNM

VBM_016_E_03/2026 Subject to technical changes and errors.